






POWERFLOW SINGLE USE SUCTION IRRIGATION INSTRUCTIONS FOR USE

Including sterile packaged PowerFlow Single-Use Suction Irrigation Devices

 **BEFORE USING THIS PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**

IMPORTANT!

1. This booklet is designed to assist using this product. It is not a reference to surgical techniques.
2. This device was designed, tested and manufactured for single use patient use only.
3. Reuse of this device may lead to its failure and subsequent patient injury.
4. Reprocessing and/or reesterilization of this device may create risk of contamination and patient infection.
5. DO NOT reuse, reprocess or reesterilize this device.

DESCRIPTION

The GENICON PowerFlow Single-Use Suction Irrigation device consists of a hand-piece equipped with two trumpet style valves, a probe with vents to produce suction and irrigation, and connecting lines of tubing, one set designed to attach to a supply of irrigation fluid, and the other designed to attach to an aspiration pump. The valves allow controlled irrigation and aspiration during a surgical procedure. The motor driven impeller pump is powered by one 9V alkaline battery (ANSI-1604A, IEC-6LR61). The battery, motor, and impeller are all located inside the hand-piece.

PROBES/ELECTROCAUTERY COMPATIBILITIES

The GENICON PowerFlow Single-Use Suction Irrigation devices are compatible with GENICON electrocautery and Stryker probes.

INDICATIONS FOR USE

The GENICON PowerFlow Single-Use Suction Irrigation device is available with an array of probe designs to facilitate lavage during laparoscopic surgery. This device has applications in laparoscopic, gynecologic, general, thoracic and urology procedures to provide suction and irrigation functions to help flush blood and tissue debris from the operative site during laparoscopy to aid visualization.

CONTRAINDICATIONS

1. The GENICON PowerFlow Single-Use Suction Irrigation device is NOT intended for contraceptive coagulation of fallopian tissue.
2. This device is not intended for used except as indicated.
3. Use of this device for hysteroscopy or for cavity distention is contraindicated.

WARNINGS AND PRECAUTIONS

1. The use of GENICON laparoscopic and endoscopic products is contraindicated whenever endoscopic surgical techniques are contraindicated for any reason. Contraindications relevant to individual GENICON products are noted in the specific sections.
2. Consult medical literature or country specific regulations for specific techniques, complications, and hazards prior to the procedure.
3. Endoscopic surgery should be performed only by physicians who are thoroughly trained in endoscopic techniques and failure modes, precautions, and corrective actions in the event of failure.
4. No modification of this equipment is allowed.
5. This device is NOT intended to generate high pressure irrigant streams in hysteroscopic procedures or encapsulated areas.
6. When not using suction, ensure that the trumpet valve providing suction does not remain pressed so as to minimize loss pneumoperitoneum.
7. As with any surgical procedure where extensive irrigation is utilized, unintended fluid uptake levels should be closely monitored by the attending physician.
8. Never pass irrigation fluid through the probe during coagulation.
9. When using the mono-polar probes, verify that the tip is not in contact with a conductive irrigation fluid.
10. Sufficient care and distance must be maintained during use to prevent arcing to other instruments.
11. This device is provided STERILE and is intended for SINGLE procedure ONLY. DISCARD AFTER USE. DO NOT RESTERILIZE.
12. Do not use any pre-sterilized products if the packaging sterile barrier is damaged or opened prior to use.
13. Care must be taken when using laparoscopic instrumentation to avoid damage to major vessels and other anatomic structures.
14. Establish and maintain adequate pneumoperitoneum to reduce the risk of injury to internal structures.
15. Properly position the patient and note anatomical landmarks to introduce devices without unintended damage.
16. Do not use excessive force or in a manner not consistent with normal instrumentation use.
17. A thorough understanding of the principles and techniques involved in laparoscopic electrosurgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel, and damage to the device or other surgical instruments.
18. When using electrosurgery tips or electrode inserts, follow all instructions, warnings, precautions, and contraindications provided with these devices.
19. Carry out electrosurgical procedures using only non-combustible gases (CO₂).
20. Verify that the devices are compatible with other products that will be used in surgery prior to the procedure.
21. Transport and store all products safely at ambient temperatures of 5°C to 30°C and relative humidity of 35% to 65% (non-condensing). Handle with care.
22. GENICON Suction Irrigation devices are DEHP / Latex free and sterilized using gamma sterilization, except as noted.
23. Dispose of all used or damaged products using standard hospital practices for biohazard control. The battery, which powers the product, may require specific/additional disposal protocols or controls. Please refer to your specific national and/or regional guidelines concerning alkaline battery disposal.
24. Do not dispose of the battery in fire, replace, recharge, or put in backwards – may explode or leak and cause personal injury.
25. Prolonged operation without fluid in the system may damage the pump motor or seal and must be avoided. Caution should be used when initially priming the system and during bag changes that the fluid is moving through the lines.
26. Do not occlude suction tip side vent holes. If side holes become occluded, negative pressure at the tip may increase and cause tissue to be pulled into the tip potentially damaged tissue.
27. The essential performance of the GENICON PowerFlow Single-Use Suction Irrigation device complies with IEC 60601-1.
28. Device performance may be compromised at suction vacuum settings greater than 50 cm Hg.

POTENTIAL COMPLICATIONS FOR ALL GENICON® PRODUCTS

Potential complications associated with the use of laparoscopic devices include but are not limited to organ or vessel damage/perforation, hemorrhage, hematoma, trocar site hernia, and septicemia/infection.

INSTRUCTIONS FOR USE

1. Inspect all components carefully for any damage that may have occurred during shipment. Inspect packaging marked as sterile, and verify that sterile barriers are intact. Do not use if sterile package is damaged or opened.
2. Open the package by using standard sterile technique.
3. Check the device for damage – do not use if the unit is defective. **Refer to picture 1**
 - A. Probe (if included)
 - B. Handpiece
 - C. Tubing
 - D. Suction Connector
 - E. Bag Spike(s) with Protectors(s)
4. To attach/tighten the probe, twist in a clockwise motion. To remove the probe, twist the probe nut in a counter clockwise motion. If probe is equipped with a mono-polar post, tighten until post points upward. **Refer to picture 2**

5. To hook up the irrigation tubing to the irrigation bag: remove the tip protector from the bag spike(s) (**Refer to picture 3, step 1**) and insert the bag spike(s) to the irrigation fluid bag (**Refer to picture 3, step 2**). If the system has two bag spikes, both should be attached to eliminate the loss of irrigation fluid.
6. Press the irrigation (Blue) trumpet valve to gravity prime the unit until fluid exits probe tip and all air has been purged. Priming must occur outside of the patient. Re-priming the unit after spiking a new saline bag may be required.
7. Attach the Suction Tubing to the Suction Canister or wall unit by inserting the blue connector located at the end of the tubing into the suction canister. **Refer to picture 4**
8. Insert shaft of probe through a laparoscopic cannula advancing gently to penetrate the cannula valve. Under appropriate visualization, advance until tip is at desired location.
9. To irrigate, press the blue colored trumpet valve. **Refer to picture 5**
10. To activate suction, press the green colored trumpet valve. **Refer to picture 6**
11. When using electrocautery probe with the GENICON PowerFlow Single-Use Suction Irrigation device attach the manufacturer's electrocautery cable to the cautery adapter provided on the back of the GENICON electrocautery probe and follow the manufacturer's instructions for use and setup.
12. **Battery Removal:** To remove the battery push the tab located on the front of the battery casing inward and pull up on the battery casing. **Refer to picture 7**

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This product may be covered by one or more U.S. patents and their foreign counter parts. Other Patents may be pending.
GENICON and the GENICON logo design are U.S. registered trademarks of GENICON Rights Reserved.

Disclaimer of Warranties on Re-Processing or Re-Sterilizing of Vendor's Disposable (Single-Patient) Products: Any reprocessing or re-sterilizing of disposable single patient devices voids Vendor's warranty. Vendor assumes no liability for reprocessed or re-sterilized products. Vendor warrants that reasonable care has been used in the manufacture of these disposable single patient devices and that they are free from defects in workmanship or materials at the time of shipment from Vendor. Vendor's sole obligation shall therefore be to repair or replace any device which it determines was defective at the time of shipment. Because no product is completely effective under all circumstances, and because the actual use and handling of these devices are beyond Vendor's control, Vendor cannot warrant for a good effect or against a bad effect in the application and use of these disposable single patient devices. The Purchaser therefore assumes all liability resulting from re-sterilization or re-processing of these disposable single patient products. Vendor therefore gives no warranty of merchantability or fitness for a particular purpose. Vendor will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to a disposable product that has been re-sterilized or re-processed. Vendor shall not be liable for incidental or consequential loss, damage or expense resulting from the use or application of these said products. This warranty is in lieu of all other warranties, whether implied, express, oral or written, and no individual or entity has the authority to vary the terms of this warranty.

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