IMPORTANT!

This booklet is designed to assist in using this product. It is not a reference for surgical techniques.

DESCRIPTION

The GENICON Reusable Clip Appliers are designed to be used on single clips that have to be loaded individually. The clip applier is designed for use during an open procedure.

INDICATIONS

The GENICON® reusable clip applier can be utilized in open procedures for occlusion of vessels and other small tissue structures, as well as for radiographic markings.

CONTRAINDICATIONS

1. Do not use the GENICON® Clip Applier where metal ligating clips would not normally be used.
2. This device is not intended for contraceptive tubal occlusion.

WARNINGS AND PRECAUTIONS

1. The GENICON® reusable clip applier must be carefully inspected prior to and after each use; because any damage may affect the safe operation of the device. Devices with moving parts should be evaluated for proper function prior to application. If any damage is suspected, contact your local representative for replacement or repair by the manufacturer.
2. The reusable clip applier is provided nonsterile and must be cleaned and sterilized prior to each use.
3. The clip applier and clips should not get in contact with electrosurgical electrodes.
4. The instrument is not intended for use as a grasper or dissector.
5. Appliers should only be used with GENICON clips
6. To prevent premature closure or release of the clip:
 - Do not squeeze applier handles before jaws are placed around vessel or tissue.
 - Ensure appliers are fully closed (handles are squeezed as close together as possible) when applying clip.
 - Inspect ligation site to ensure secure application and hemostasis.
7. To prevent damage to applier jaws:
 - Avoid introducing applier tips at an angle into the cartridge rack, this may misalign jaws.
 - Inspect applier and jaws before each use.
 - Use only GENICON branded clips.

Attention – Federal (USA) law restricts this device to sale by or on the order of a physician.

Attention – Clip applier will arrive in open package. Inspect product thoroughly if package is damaged.

Attention – Product must be cleaned and sterilized PRIOR to EACH use. Please see cleaning instructions below.

Important – This booklet is designed to provide instructions for the use of an open clip applier. It is not a reference for surgical techniques. This product should only be used by a properly trained surgeon. Before using this product, read the following information thoroughly. Failure to properly follow the instructions may lead to serious surgical injuries.

Important:

- This clip is not intended for Fallopian Tube Ligation.
- As a non-magnetic marking clip, it has low magnetic susceptibility and is safe for use in the MRI environment. Positional errors and artifacts may increase dramatically with increasing static magnetic field strength and/or increasing field of view size. Artifact size may vary greatly for different MRI pulse sequences.
- The titanium clip is biologically inert and can remain in body without consequences.

Contents

DESCRIPTION	SIZE	PRODUCT CODE	CLIP SIZE
Straight Jaw	14.6 cm	330-005-000	Small
Straight Jaw	26.7 cm	330-005-001	Small
Straight Jaw	19.1 cm	330-005-002	Medium
Straight Jaw	26.7 cm	330-005-003	Medium
Angled Jaw	26.7 cm	330-005-004	Medium/Large
Angled Jaw	26.7 cm	330-005-005	Large

Instructions

1. Inspect the integrity of the individual protection (packaging) before initial use. The package is not reusable – discard of it appropriately after shipment.
2. Open the packaging. *product DOES NOT come sterilized. All applicators MUST be CLEANED AND STERILIZED PRIOR to EACH use*
3. Place clip cartridge firmly on a stable surface and hold clip cartridge in place while removing clips.
4. Select applicator according to clip size.
5. Introduce the jaws of the applicator vertically into the cartridge clip rack.
6. Push the applicator all the way into the cartridge rack and withdraw to remove the clip.
7. Place the loaded clip and jaws of the applicator around the vessel or tissue to be occluded.
8. Firmly squeeze the handles together closing the clip. The clip closes tip to tip preventing vessel slippage before fully closing.
9. Release the handle and remove applicator. If handle is not fully decompressed clip may remain on applicator resulting in improper ligation, injury or other serious surgical complications.

Cleaning instructions

The user must ensure that cleaning and sterilization are conducted in accordance with the appropriate guidelines, standards or national Health Authorities' requirements.

Sterilization of the device may be accomplished by steam. When carrying out the steam sterilization – the sterilization has to be done according to a validated steam sterilization procedure (e.g. in a sterilizer according to EN ISO 11137-1 (2006) EN ISO 11137-2 (2006) and EN ISO 11137-3 (2006)).

GENICON does not recommend the device be sterilized by “Flash” or chemical sterilization.

Applicator may also be placed within a rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

Minimal recommended sterilization parameters are as follows:

Sterilization Method	Temp	Minimum Exposure Time	
		Wrapped	Sterile Container
Pre-Vacuum	270 – 275 °F	4 minutes	4 minutes

Storage

Store the clip applicator in a cool dry area, with good ventilation away from areas of extreme temperature and humidity..

Shelf Life

Do not use if package is damaged.

The manufacturer declines responsibility for problems resulting from improper use and/or improper cleaning/sterilization.

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REMAINDER OF PAGE INTENTIONALLY
LEFT BLANK

This product may be covered by one or more U.S. patents and their foreign counter parts. Other Patents may be pending.
GENICON and the GENICON logo design are U.S. registered trademarks of GENICON Rights Reserved.

Disclaimer of Warranties on Re-Processing or Re-Sterilizing of Vendor's Disposable (Single-Patient) Products: Any reprocessing or re-sterilizing of disposable single patient devices voids Vendor's warranty. Vendor assumes no liability for reprocessed or re-sterilized products. Vendor warrants that reasonable care has been used in the manufacture of these disposable single patient devices and that they are free from defects in workmanship or materials at the time of shipment from Vendor. Vendor's sole obligation shall therefore be to repair or replace any device which it determines was defective at the time of shipment. Because no product is completely effective under all circumstances, and because the actual use and handling of these devices are beyond Vendor's control, Vendor cannot warrant for a good effect or against a bad effect in the application and use of these disposable single patient devices. The Purchaser therefore assumes all liability resulting from re-sterilization or re-processing of these disposable single patient products. Vendor therefore gives no warranty of merchantability or fitness for a particular purpose. Vendor will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to a disposable product that has been re-sterilized or re-processed. Vendor shall not be liable for incidental or consequential loss, damage or expense resulting from the use or application of these said products. This warranty is in lieu of all other warranties, whether implied, express, oral or written, and no individual or entity has the authority to vary the terms of this warranty.

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**REMAINDER OF PAGE INTENTIONALLY
LEFT BLANK**

	EN ES IT	Manufacturer Fabricante Produttore		EN ES IT	Serial Number Numero de serie Numero di serie
	EN ES IT	Latex Free Sin látex Privo di lattice		EN ES IT	Do not Use if Package is Damaged No usar si el envoltorio está dañado Non utilizzare se l'imballaggio è danneggiato
	EN ES IT	Manufacture Date Fecha de fabricacion Data di produzione		EN ES IT	Authorized Representative in EU Representante autorizado en la Unión Europea Rappresentante autorizzato UE
	EN ES IT	Caution Cuidado Cautela		EN ES IT	Authorized For Sale or Use by Physician Only Autorizado para la venta o para el uso sólo por parte de un medico Autorizzato solo per vendita o uso da personale medico
INTENTIONALLY LEFT BLANK				EN ES IT	See Instructions for Use Consultar las instrucciones de uso Vedere le istruzioni per l'uso

	6869 Stapoint Court Suite #114 Winter Park, Florida 32792 USA TEL: 1-407-657-4851 FAX: 1-407-677-9773 www.geniconendo.com info@geniconendo.com		European Authorized Representative Genicon Europe 44 Airedale Avenue Chiswick, London United Kingdom W4 2NW TEL: 44-203-286-3414
---	---	--	---