






SINGLE USE PNEUMO-NEEDLE INSTRUCTIONS FOR USE

Including sterile packaged Single-Use Pnuemo-Needle

BEFORE USING THIS PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

IMPORTANT!

1. This booklet is designed to assist using this product. It is not a reference to surgical techniques.
2. This device was designed, tested and manufactured for single patient use only.
3. Reuse of this device may lead to its failure and subsequent patient injury.
4. Reprocessing and/or resterilization of this device may create risk of contamination and patient infection.
5. DO NOT reuse, reprocess or resterilize this device

DESCRIPTION

The GENICON Pnuemo-Needle is a sterile, disposable Veress needle which is available in 120mm and 150mm length. The device is equipped with a spring-loaded, round-tipped obturator. In addition, there is a “slide switch” which permits easy ON-OFF control of the gas flow. The most proximal end contains a male luer lock connector for secure CO₂ gas line connection.

INDICATIONS FOR USE

The GENICON Pnuemo-Needle is intended for percutaneous insertion into the peritoneal cavity for the purpose of insufflation with carbon dioxide to establish pneumoperitoneum prior to the placement of trocars during laparoscopic procedures.

CONTRAINDICATIONS

- The use of GENICON laparoscopic and endoscopic products is contraindicated whenever endoscopic surgical techniques are contraindicated for any reason. Contraindications relevant to individual GENICON products are noted in the specific sections.

WARNINGS AND PRECAUTIONS

- Endoscopic surgery should be performed only by physicians who are thoroughly trained in endoscopic techniques and failure modes, precautions, and corrective actions in the event of failure.
- GENICON disposable products are intended for single patient use only – DO NOT RESTERILIZE.
- Do not use any pre-sterilized products if the packaging sterile barrier is damaged
- Consult medical literature or country specific regulations for specific techniques, complications, and hazards prior to procedure.
- Care must be taken when using laparoscopic instrumentation to avoid damage to major vessels and other anatomic structures.
- Establish and maintain adequate pneumoperitoneum to reduce the risk of injury to internal structures.
- Properly position the patient and note anatomical landmarks to introduced devices without unintended damage.
- Do not use excessive force or in a manner not consistent with normal instrumentation use.
- A thorough understanding of the principles and techniques involved in laparoscopic electrosurgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel, and damage to the device or other surgical instruments.
- Carry out electrosurgical procedures using only non-combustible gases (CO₂).
- Verify that the devices are compatible with other products that will be used in surgery prior to the procedure.
- Store all products safely in a climate controlled environment and handle with care.
- All GENICON devices are latex free and sterilized using gamma sterilization, except as noted.
- Dispose of all used or damaged products using standard hospital practices for biohazard control.

POTENTIAL COMPLICATIONS FOR ALL GENICON® PRODUCTS

Potential complications associated with the use of laparoscopic devices include but are not limited to organ or vessel damage/perforation, hemorrhage, hematoma, trocar site hernia, and septicemia/infection.

INSTRUCTIONS FOR USE

1. Inspect sterile packaging, and verify that sterile barriers are intact. Do not use if package is compromised.
2. Remove needle protector. **Refer to picture 1**
3. Inspect all components carefully for damage during shipment, handling or after maintenance. **Refer to picture 2**
 - A. CO₂ Luer Port/Valve
 - B. Indicator
 - C. Grip
 - D. Outer Needle
 - E. Inner Blunt Tip

4. The insufflation (CO₂) valve is packaged in the CLOSED position (red slide valve fully depressed). To insufflate the abdomen, attach the CO₂ supply to the luer connection, and press the green colored slide valve fully. **Refer to picture 3** To stop insufflation and close the valve, press the red colored slide valve fully. **Refer to picture 4**
5. Using laparoscopic technique, select entry site advance the needle through the skin and fascia.
6. Attach a syringe to the Luer fitting. **Refer to picture 5**
7. Open the airflow switch, and draw back the syringe plunger to aspirate through the needle. **Refer to picture 6**
Presence of blood or other liquid indicates potential perforation of organs or vessels.
8. If there is NO aspiration, the insufflation line may be connected to obtain pneumoperitoneum. **Refer to picture 7**
9. Once sufficient insufflation is obtained, position the airflow switch to the "OFF" position and remove the needle when desired.

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This product may be covered by one or more U.S. patents and their foreign counter parts.
Other Patents may be pending.

GENICON and the GENICON logo design are U.S. registered trademarks of GENICON
Rights Reserved.

Disclaimer of Warranties on Re-Processing or Re-Sterilizing of Vendor's Disposable (Single-Patient) Products: Any reprocessing or re-sterilizing of disposable single patient devices voids Vendor's warranty. Vendor assumes no liability for reprocessed or re-sterilized products. Vendor warrants that reasonable care has been used in the manufacture of these disposable single patient devices and that they are free from defects in workmanship or materials at the time of shipment from Vendor. Vendor's sole obligation shall therefore be to repair or replace any device which it determines was defective at the time of shipment. Because no product is completely effective under all circumstances, and because the actual use and handling of these devices are beyond Vendor's control, Vendor cannot warrant for a good effect or against a bad effect in the application and use of these disposable single patient devices. The Purchaser therefore assumes all liability resulting from re-sterilization or re-processing of these disposable single patient products. Vendor therefore gives no warranty of merchantability or fitness for a particular purpose. Vendor will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to a disposable product that has been re-sterilized or re-processed. Vendor shall not be liable for incidental or consequential loss, damage or expense resulting from the use or application of these said products. This warranty is in lieu of all other warranties, whether implied, express, oral or written, and no individual or entity has the authority to vary the terms of this warranty.

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