



Simply perfect from  
every point of view.



## Bladeless Visual Tip Trocar and Cannula System

BY THE SOCIETY OF  
LAPAROENDOSCOPIC SURGEONS

New Ideas for Modern Health Care

## Access Ports: Bladeless Visual Tip Trocar and Cannula System

The distinctively engineered visual tip greatly reduces wound defect size as well as insertion force into the abdominal cavity.

- ✓ Bladeless tip dissects without cutting and allows for visualization during insertion
- ✓ Injuries due to blind entry and a sharp blade are virtually eliminated with the use of the GeniLook<sup>™</sup> Trocar System
- ✓ Universal handle can be used as a "pistol grip" or like a traditional Trocar handle when using a scope
- ✓ Cannula incorporates a beveled and angled distal tip
- ✓ Dual seal valve technology minimizes loss of insufflation pressure, maintaining pneumoperitoneum



High-visualization bladeless tip allows for visualization during insertion

### One Cannula

### Two Cannula

Diameter	Length	Qty/Box	Threaded	Threaded
5mm	70mm	5	200-005-051	200-005-052
5mm	100mm	5	200-005-151	200-005-152
5mm	150mm	5	200-005-251	200-005-252
10mm	100mm	5	200-010-151	200-010-152
10mm	150mm	5	200-010-251	---
12mm	100mm	5	200-012-151	---
12mm	150mm	5	200-012-251	---
15mm	100mm	5	200-015-151	---
15mm	150mm	5	200-015-251	---

Every 15mm Cannula comes with a 5mm Vortex Reducer Valve