



When distance  
is never an issue.



**GENISurge**<sup>®</sup>



BY THE SOCIETY OF  
LAPAROENDOSCOPIC SURGEONS

Disposable Monopolar

New Ideas for Modern Health Care



Curved Metzenbaum Scissors



Straight Scissors



Fully Insulated Scissors  
Full length insulation protects tissue from electrocautery arcs



Curved Maryland Dissector



Atraumatic Fenestrated Grasper




Babcock Grasper



Maxi Grip Fenestrated Grasper

**GeniSurge 33cm Monopolar instruments offer a color coded knob-to-jaw system making it easier to differentiate inside the operating room.**



Scissors



Grasper



Maxi Grip



Dissector



Babcock

## Disposable Monopolar

Tip Description	Qty/Box	Shaft Length	Ratcheting	Non-Ratcheting
Curved Maryland Dissector	10	33cm	533-005-950	533-005-905
Atraumatic Fenestrated Grasper	10	33cm	533-005-920	533-005-915
Babcock Grasper	10	33cm	533-005-925	533-005-930
Maxi Grip Fenestrated Grasper	10	33cm	533-005-945	533-005-940
Curved Metzenbaum Scissors	10	33cm		533-005-910
Straight Scissors	10	33cm		533-005-935
Fully Insulated Scissors	10	33cm		533-105-910