



## Sleek Vision all around



## Natura Optical Bladeless Trocar

## Access Ports: Optical Bladeless Trocar and Cannula System

The distinctively engineered visual tip greatly reduces wound defect size as well as insertion force into the abdominal cavity.

High-visualization bladeless tip allows for visualization during insertion



- ✓ Dual seal valve technology minimizes loss of insufflation pressure, maintaining pneumoperitoneum
- ✓ Injuries due to blind entry and a sharp blade are virtually eliminated
- ✓ Ergonomic wing handle provides superior comfort and stability
- ✓ Cannula incorporates a beveled and angled distal tip

Diameter	Length	Qty/Box	One Cannula	Two Cannula
			Threaded	Threaded
5mm	70mm	5	210-005-051	210-005-052
5mm	100mm	5	210-005-151	210-005-152
10mm	100mm	5	210-010-151	210-010-152
12mm	100mm	5	210-012-151	---
15mm	100mm	5	210-015-151	---

Configurations with Smooth Cannula available for order in minimum quantities  
Every 15mm Cannula comes with a 5mm Vortex Reducer Valve