



Xceed Your Expectations



Instrumentation

New Ideas for Modern Health Care



Compatible with Surgical Innovation®s
101-41000, 101-42000, 101-43000,
101-43500, and 101-44000 handles

Description	Qty/Box	Size	Product Code
Mini Curved Metz. Scissor Insert/Shaft	10	5mm x 33cm	590-005-005
Standard Curved Metz. Scissor Insert/Shaft	10	5mm x 33cm	590-005-010
Standard Curved Metz. Scissor Insert/Shaft	10	5mm x 45cm	590-005-013

Hand Sharpened

Honed Blades provide superior cutting through the entire stroke

Fitted Surgical Stainless Steel Blades

Each blade is carefully balanced and matched to a precise cutting arc to maximize visualization and functionality of the blade interaction

Enhanced Cutting

at the distal tip maintains accurate and precise tip approximation for improved micro dissection



Compatible with Microline® 3904, 3905, 3908, 3914, 3915, 3918, 3924, 3925, 3934, and 3935 handles

Description	Qty/Box	Size	Product Code
Curved Metzenbaum Scissor Tip	10	5mm x 5.6cm	590-001-001

The GENICON X-Surge Tips and Shafts are intended for use with electrosurgical handles in compliance with safety standards. IEC 60601-1, IEC 60601-1-2, and IEC 60601-2-2 including applicable national/regional differences. Follow all safety precautions.

Compatibility of parts from different manufacturers must be verified prior to use. SURGICAL INNOVATIONS® is a registered trademark of SURGICAL INNOVATIONS®, LLC.

MICROLINE® is a registered trademark of MICROLINE SURGICAL®, Inc.

XSURGE™ is a trademark of GENICON®, Inc. XSURGE™ is neither associated with nor sponsored by SURGICAL INNOVATIONS®, LLC and MICROLINE SURGICAL®, Inc.

6869 Stapoint Court #114 | Winter Park, Florida USA 32792 | Tel: 1.407.657.4851 Fax: 1.407.677.9773 | www.geniconendo.com

© 1998-2017 GENICON® All rights Reserved. US and foreign patents issued or pending. Visit www.geniconendo.com for a list of trademarks and patents